

Product datasheet for **AM08096LE-N**

Rat IgM Mu chain specific F(ab')₂ Fragment Mouse Monoclonal Antibody [Clone ID: M 2A1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgM Mu chain specific F(ab') ₂ Fragment Mouse Monoclonal Antibody [Clone ID: M 2A1]
Clone Name:	M 2A1
Applications:	FN
Recommended Dilution:	B Cell Stimulation: ≤5 µg per mL of tissue culture media.
Reactivity:	Rat
Host:	Mouse
Immunogen:	Rat IgM Hybridoma.
Isotype:	IgG1
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)