

Product datasheet for **AM08099BT-N**

Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]
Clone Name:	SB87a
Applications:	ELISA, FC
Recommended Dilution:	ELISA: 1/5,000-1/20,000. Flow Cytometry.
Reactivity:	Rabbit
Host:	Rat
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)