

Product datasheet for **AM08099HR-N**

Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]
Clone Name:	SB87a
Applications:	ELISA
Recommended Dilution:	ELISA: 1/4,000-1/8,000.
Reactivity:	Rabbit
Host:	Rat
Isotype:	IgG2a
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)