

## Product datasheet for **AM08099PU-N**

### Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit Ig (IgG + IgM) Rat Monoclonal Antibody [Clone ID: SB87a]
Clone Name:	SB87a
Applications:	ELISA, FC
Recommended Dilution:	<b>ELISA:</b> 1/4,000-1/8,000. <b>Flow Cytometry.</b>
Reactivity:	Rabbit
Host:	Rat
Immunogen:	Rabbit IgG and IgM
Isotype:	IgG2a
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)