

Product datasheet for **AM08102AP-N**

Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9] |
| Clone Name: | 2A9 |
| Applications: | ELISA |
| Recommended Dilution: | ELISA: 1/2,000-1/4,000. |
| Reactivity: | Rabbit |
| Host: | Mouse |
| Isotype: | IgG1 |
| Formulation: | 50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0 containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase |
| Conjugation: | AP |
| Storage: | Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |



[View online »](#)