

## Product datasheet for **AM08102BT-N**

### Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]
Clone Name:	2A9
Applications:	ELISA, FC
Recommended Dilution:	<b>ELISA:</b> 1/5,000-1/20,000. <b>Flow Cytometry.</b>
Reactivity:	Rabbit
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)