

Product datasheet for AM08102FC-N

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Product data:

Product Type: Secondary Antibodies

Product Name: Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Clone Name: 2A9
Applications: FC

Recommended Dilution: Flow Cytometry: $< / = 3 \mu g/10e6$ cells.

Reactivity: Rabbit
Host: Mouse
Isotype: IgG1

Formulation: PBS containing 0.09% Sodium Azide as preservative.

Label: FITC

State: Liquid purified Ig fraction.

Label: Fluorescein Isothiocyanate Isomer 1

Concentration: lot specific

Conjugation: FITC

Storage: Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.