

## Product datasheet for AM08102HR-N

## 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436

OriGene Technologies, Inc.

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

## Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

**Product data:** 

**Product Type:** Secondary Antibodies

**Product Name:** Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Clone Name: 2A9

**Applications:** ELISA

Recommended Dilution: ELISA: 1/4,000-1/8,000.

Reactivity: Rabbit Host: Mouse

Immunogen: Recombinant Rabbit IgG Fc Fragment

Isotype: IgG1

**Formulation:** 50% Glycerol/50% PBS, pH 7.4 without preservatives.

Label: HRP

State: Liquid purified Ig fraction. Label: Horseradish Peroxidase

Conjugation: HRP

**Storage:** Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.