

Product datasheet for **AM08102PU-N**

Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]
Clone Name:	2A9
Applications:	ELISA, FC
Recommended Dilution:	ELISA: 1/4,000-1/8,000. Flow Cytometry. FLISA.
Reactivity:	Rabbit
Host:	Mouse
Isotype:	IgG1
Formulation:	Borate Buffered Saline, pH 8.2 State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)