

Product datasheet for **AM08102RP-N**

Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rabbit IgG Mouse Monoclonal Antibody [Clone ID: 2A9]
Clone Name:	2A9
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 0.3 µg/10e6 cells.
Reactivity:	Rabbit
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)