

## Product datasheet for **AM08103RP-N**

### CD1.1 Mouse Monoclonal Antibody [Clone ID: 76-7-4]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	76-7-4
Applications:	FC
Recommended Dilution:	<b>Flow Cytometry:</b> < / = 0.2 µg/10e6 cells. (Ref.1)
Reactivity:	Porcine
Host:	Mouse
Isotype:	IgG2a
Clonality:	Monoclonal
Specificity:	This antibody is specific to Porcine CD1 (Mr 43-49 kDa).
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">Entrez Gene 396785 Pig</a> <a href="#">Q9XS72</a>
Background:	All CD1 molecules, except CD1e, are cell surface glycoproteins that are structurally related to the MHC molecules, however, in distinction, CD1 proteins are essentially non polymorphic. CD1 has considerable structural homology with both MHC class I and class II molecules, and CD1 molecules are involved in T cell activation. In contrast to MHC, however, CD1 molecules appear to present predominantly non peptide molecules originating from lipids and glycolipids.



[View online »](#)

Synonyms: T6/Leu-6, hTa1