

Product datasheet for **AM08125BT-N**

Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-4]
Clone Name:	M-4
Applications:	FC, IF
Recommended Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells. Indirect Immunofluorescence staining. (Ref.1)
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgM
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)