

Product datasheet for **AM08125PU-N**

Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-4]
Clone Name:	M-4
Applications:	FC, FN, IF
Recommended Dilution:	Flow Cytometry: $\leq 1 \mu\text{g}/10^6$ cells. Indirect Immunofluorescence staining. (Ref.1) <i>In ovo</i> suppression of Chicken B lymphocytes. <i>In vitro</i> activation of Chicken peripheral B lymphocytes.
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgM
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)