

Product datasheet for **AM08126FC-N**

Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]


Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]
Clone Name:	M-1
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 1 µg/10e6 cells. (Ref.1)
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)

Product images:



Immunofluorescent Staining: Chicken peripheral blood mononuclear cells were double-stained with Mouse anti-Chicken IgM-FITC and Mouse anti-Chicken CD3-PE. Small lymphocytes were then gated and analyzed on a FACScan(TM) flow cytometer (BDIS, San Jose, CA). Amount Used: 1 µg/10e6 cells.