

## Product datasheet for **AM08126PU-N**

### Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]
Clone Name:	M-1
Applications:	FC, IP
Recommended Dilution:	<b>Flow Cytometry.</b> <b>Immunoprecipitation.</b> (Ref.1)
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG2b
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)