

Product datasheet for **AM08126RP-N**

Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgM (Mu heavy chain specific) Mouse Monoclonal Antibody [Clone ID: M-1]
Clone Name:	M-1
Applications:	FC
Recommended Dilution:	Flow Cytometry: < / = 0.2 µg/10e6 cells. (Ref.1)
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)