

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## Product datasheet for AM08127FC-N

## Chicken IgG / Chicken IgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1]

## **Product data:**

| Product Type:         | Secondary Antibodies  |
|-----------------------|---|
| Product Name:         | Chicken lgG / Chicken lgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1]  |
| Clone Name:           | G-1   |
| Applications:         | FC, IF  |
| Recommended Dilution: | Flow Cytometry: < / = 1 µg/10e6 cells.<br>Immunofluorescence.   |
| Reactivity:           | Chicken   |
| Host:                 | Mouse   |
| lsotype:              | lgG1  |
| Formulation:          | PBS containing 0.09% Sodium Azide as preservative.<br>Label: FITC<br>State: Liquid purified Ig fraction.<br>Label: Fluorescein Isothiocyanate Isomer 1  |
| Concentration:        | lot specific  |
| Conjugation:          | FITC  |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.<br>This product is photosensitive and should be protected from light.<br>Avoid repeated freezing and thawing. |



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US