

Product datasheet for **AM08127PU-N**

Chicken IgG / Chicken IgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgG / Chicken IgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1]
Clone Name:	G-1
Applications:	ELISA, FC, IF, IHC
Recommended Dilution:	ELISA. Flow Cytometry. Immunofluorescence. Immunohistochemistry on Frozen Sections.
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)