

Product datasheet for **AM08127RP-N**

Chicken IgG / Chicken IgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1]

Product data:

| | |
|-----------------------|---|
| Product Type: | Secondary Antibodies |
| Product Name: | Chicken IgG / Chicken IgY (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: G-1] |
| Clone Name: | G-1 |
| Applications: | FC, IF |
| Recommended Dilution: | Flow Cytometry: < / = 0.2 µg/10e6 cells. |
| Reactivity: | Chicken |
| Host: | Mouse |
| Isotype: | IgG1 |
| Formulation: | PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin |
| Concentration: | lot specific |
| Conjugation: | PE |
| Storage: | Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. |



[View online »](#)