

## Product datasheet for **AM08128BT-N**

### Chicken IgA (alpha chain specific) Mouse Monoclonal Antibody [Clone ID: A-1]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Chicken IgA (alpha chain specific) Mouse Monoclonal Antibody [Clone ID: A-1]
Clone Name:	A-1
Applications:	ELISA, FC, IF
Recommended Dilution:	<b>ELISA.</b> <b>Flow Cytometry:</b> < / = 1 µg/10e6 cells. <b>Immunofluorescence.</b>
Reactivity:	Chicken
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)