

Product datasheet for **AM08145BT-N**

Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]
Clone Name:	SA-DA4
Applications:	ELISA, IF
Recommended Dilution:	Indirect Immunofluorescent staining of IgM+ human B lymphocytes: $< / = 1 \mu\text{g}/10\text{e}6$ cells. ELISA: 1/5,000-1/10,000.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)