

Product datasheet for **AM08145PU-N**

Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]
Clone Name:	SA-DA4
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Indirect Immunofluorescent staining of IgM+ human B lymphocytes: $\leq 1 \mu\text{g}/10^6$ cells. ELISA. Immunohistochemistry on Frozen Sections. Western blotting. Dot- and Slot-Immunoblotting. Immunocytochemistry.
Reactivity:	Human
Host:	Mouse
Immunogen:	Human IgM myeloma protein.
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)