

## Product datasheet for **AM08145RP-N**

### Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: SA-DA4]
Clone Name:	SA-DA4
Applications:	IF
Recommended Dilution:	<b>Indirect Immunofluorescent staining</b> of IgM+ human B lymphocytes: $< / = 0.2 \mu\text{g}/10\text{e}6$ cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)