

Product datasheet for **AM08146LE-N**

Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: UHB]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain specific) Mouse Monoclonal Antibody [Clone ID: UHB]
Clone Name:	UHB
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Indirect Immunofluorescent staining of IgM+ human B lymphocytes. ELISA. Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry on Frozen Sections. Immunocytochemistry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG3
Formulation:	PBS without preservatives. State: Low Endotoxin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.



[View online »](#)