

Product datasheet for **AM08148BT-N**

Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	ELISA, IF
Recommended Dilution:	Indirect Immunofluorescent staining of IgD+ human B lymphocytes: $\leq 1 \mu\text{g}/10^6$ cells. ELISA: 1/5,000-1/10,000.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)