

Product datasheet for **AM08148FC-N**

Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	IF
Recommended Dilution:	Indirect Immunofluorescent staining of IgD+ human B lymphocytes: < / = 1 µg/10e6 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)