

## Product datasheet for **AM08148HR-N**

### Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/1,000-1/2,000.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)