

Product datasheet for **AM08148PU-N**

Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Indirect Immunofluorescent staining of IgD+ human B lymphocytes: < / = 1 µg/10e6 cells. ELISA. 1 µg/mL Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry on Frozen Sections. Immunocytochemistry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)