

## Product datasheet for **AM08148RP-N**

### Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgD (Delta chain specific) Mouse Monoclonal Antibody [Clone ID: IADB6]
Clone Name:	IADB6
Applications:	IF
Recommended Dilution:	<b>Indirect Immunofluorescent staining</b> of IgD+ human B lymphocytes: < / = 0.2 µg/10e6 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2a
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction. Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)