

Product datasheet for **AM08149AP-N**

Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: JDC-10]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: JDC-10]
Clone Name:	JDC-10
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000-1/4,000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Purified human IgG1 myeloma protein
Isotype:	IgG1
Formulation:	50mM Tris/1mM MgCl ₂ /50% Glycerol, pH 8.0, containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)