

Product datasheet for **AM08149BT-N**

Human IgG (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: JDC-10]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (gamma chain specific) Mouse Monoclonal Antibody [Clone ID: JDC-10]
Clone Name:	JDC-10
Applications:	ELISA, IF
Recommended Dilution:	Indirect Immunofluorescent staining of IgG+ human B lymphocytes: < / = 1 µg/10e6 cells. ELISA: 1/5,000-1/10,000. Flow Cytometry: 1 µg / 10e6 cells Western Blotting
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)