

## Product datasheet for **AM08150AC-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]
Clone Name:	H2
Applications:	IF
Recommended Dilution:	<b>Indirect Immunofluorescent staining</b> of IgG+ human B lymphocytes: < / = 0.1 µg/10e6 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: APC State: Liquid purified Ig fraction. Label: Allophycocyanin
Concentration:	lot specific
Conjugation:	APC
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)