

## Product datasheet for **AM08150AP-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]
Clone Name:	H2
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/2000-1/8,000.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	50mM Tris/1mM MgCl <sub>2</sub> /50% Glycerol, pH 8.0, containing 0.09% Sodium Azide as preservative. Label: AP State: Liquid purified Ig fraction. Label: Alkaline Phosphatase
Conjugation:	AP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)