

Product datasheet for **AM08150FC-N**

Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]
Clone Name:	H2
Applications:	IF
Recommended Dilution:	FLISA: 1/200-1/600. Indirect Immunofluorescent staining of IgG+ human B lymphocytes: $\leq 1 \mu\text{g}/10^6$ cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)