

Product datasheet for **AM08150HR-N**

Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]
Clone Name:	H2
Applications:	ELISA
Recommended Dilution:	ELISA: 1/4,000-1/8,000 (See also Ref.1).
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives Label: HRP State: Liquid purified Ig fraction Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)