

## Product datasheet for **AM08150PU-N**

### Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Mouse Monoclonal Antibody [Clone ID: H2]
Clone Name:	H2
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	<b>Indirect Immunofluorescent staining</b> of IgG+ human B lymphocytes: $\leq 1 \mu\text{g}/10^6$ cells. <b>ELISA.</b> <b>Western blotting.</b> <b>Dot- and Slot-Immunoblotting.</b> <b>Immunohistochemistry on Frozen Sections.</b> <b>Immunocytochemistry.</b>
Reactivity:	Human
Host:	Mouse
Isotype:	IgG2b
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)