

Product datasheet for **AM08151FC-N**

Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]
Clone Name:	4E3
Applications:	IF
Recommended Dilution:	FLISA: 1:200-1:400 Immunofluorescence < / = 1 µg/10e6 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	100 mM borate buffered saline, pH 8.2 with 0.09% sodium azide as preservative. Label: FITC State: Liquid purified Ig fraction. Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C. This product is photosensitive and should be protected from light.



[View online »](#)