

Product datasheet for **AM08151HR-N**

Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]
Clone Name:	4E3
Applications:	ELISA
Recommended Dilution:	ELISA: 1/2,000-1/8,000.
Reactivity:	Human
Host:	Mouse
Immunogen:	Unknown.
Isotype:	IgG1
Formulation:	50% Glycerol/50% PBS, pH 7.4 without preservatives. Label: HRP State: Liquid purified Ig fraction. Label: Horseradish Peroxidase
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)