

Product datasheet for **AM08151PU-N**

Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]
Clone Name:	4E3
Applications:	ELISA, IF
Recommended Dilution:	ELISA: $\leq 5 \mu\text{g/ml}$ (Ref.1) Immunofluorescence: $\leq 1 \mu\text{g}/10^6$ cells. Flow Cytometry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)