

Product datasheet for **AM08151RP-N**

Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG1 (gamma 1 chain specific) Mouse Monoclonal Antibody [Clone ID: 4E3]
Clone Name:	4E3
Applications:	ELISA, IF
Recommended Dilution:	ELISA. FLISA ($\leq 1.0 \mu\text{g/ml}$).
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent. Label: PE State: Liquid purified Ig fraction Label: R-Phycoerythrin
Concentration:	lot specific
Conjugation:	PE
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)