

## Product datasheet for **AM08152BT-N**

### Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]
Clone Name:	31-7-4
Applications:	ELISA, IF
Recommended Dilution:	<b>ELISA:</b> 1/5,000-1/20,000. <b>Immunofluorescence</b> < / = 1 µg/10e6 cells.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	PBS containing 0.09% Sodium Azide as preservative. Label: Biotin State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Biotin
Storage:	Store the antibody undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)