

Product datasheet for AM08152HR-N

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]

Product data:

Product Type: Secondary Antibodies

Product Name: Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]

Clone Name: 31-7-4
Applications: ELISA

Recommended Dilution: ELISA: 1/2,000-1/8,000.

Reactivity: Human Host: Mouse

Immunogen: Human IgG₂ myeloma protein

Isotype: IgG1

Formulation: 50% Glycerol/50% PBS, pH 7.4 without preservatives.

Label: HRP

State: Liquid purified Ig fraction. Label: Horseradish Peroxidase


Conjugation: HRP

Storage: Store undiluted at 2-8°C for one month or in (aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.



Product images:



ELISA plate was coated with purified human IgG1, IgG2, IgG3, and IgG4. Immunoglobulins were detected with serially diluted Mouse Anti-Human IgG2 Fc-HRP (AM08152HR-N).