

Product datasheet for **AM08152PU-N**

Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG2 (Fc specific) Mouse Monoclonal Antibody [Clone ID: 31-7-4]
Clone Name:	31-7-4
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	Fluorescent-linked-immunosorbent assay (FLISA) ELISA. Western blotting. Dot- and Slot-Immunoblotting. Immunohistochemistry (frozen Sections). Immunocytochemistry.
Reactivity:	Human
Host:	Mouse
Isotype:	IgG1
Formulation:	100 mM Borate Buffered Saline, pH 8.2. No preservatives or amine-containing buffer salts added. State: Purified State: Liquid purified Ig fraction.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)