

## Product datasheet for AM08154AC7-N

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## CD20 (MS4A1) Mouse Monoclonal Antibody [Clone ID: B-Ly1]

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: B-Ly1
Applications: FC

Recommended Dilution: Flow Cytometry: 10 µl/10e6 cells.

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Specificity:** This antibody recognizes CD20.

**Formulation:** PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.

Label: APC-Cy7

State: Liquid purified Ig fraction. Label: Allophycocyanin-Cyanine 7

**Conjugation:** APC-Cy7

**Storage:** Store the antibody undiluted at 2-8°C.

DO NOT FREEZE!

This product is photosensitive and should be protected from light.

Avoid repeated freezing and thawing.

Stability: Shelf life: one year from despatch.

Gene Name: membrane spanning 4-domains A1

Database Link: Entrez Gene 931 Human

P11836

**Background:** CD20 is a non glycosylated protein with a molecular weight of 35 or 37 kDa depending on the

degree of phosphorylation. Although not a member of the tetraspanin superfamily of cell surface receptors, it crosses the cell membrane four times. The CD20 antigen is present on human pre B lymphocytes and on B lymphocytes at all stages of maturation, except on plasma cells. Low level expression of the CD20 antigen has been detected on normal T lymphocytes. The CD20 molecule is involved in regulation of B cell differentiation,

presumably via its reported function as a Ca++ channel subunit. (Ref-1-4)







**Synonyms:** MS4A1, Leu-16, Bp35, B-cell marker