

## Product datasheet for **AM08206AC-N**

### Mouse IgM Isotype Control

#### Product data:

|               |  |
|---------------|--|
| Product Type: | Isotype Controls   |
| Clone Name:   | 11E10  |
| Applications: | FC   |
| Host:         | Mouse  |
| Isotype:      | IgM  |
| Clonality:    | Monoclonal   |
| Specificity:  | This antibody is specific to Lipopolysaccharide.   |
| Buffer:       | 1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.   |
| Conjugation:  | APC  |
| Storage:      | Store the antibody undiluted at 2-8°C.<br><b>DO NOT FREEZE!</b><br>This product is photosensitive and should be protected from light.<br>Avoid repeated freezing and thawing.<br>Shelf life: one year from despatch. |



[View online »](#)