

Product datasheet for **AM08206LE-N**

Mouse IgM Isotype Control

Product data:

Product Type:	Isotype Controls
Clone Name:	11E10
Applications:	E, FC
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	This antibody is specific to Lipopolysaccharide.
Buffer:	PBS without preservatives.
Concentration:	0.5 mg/ml
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)