

## Product datasheet for **AM08206PU-N**

### Mouse IgM Isotype Control

#### Product data:

Product Type:	Isotype Controls
Clone Name:	11E10
Applications:	E, FC
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	This antibody is specific to Lipopolysaccharide.
Buffer:	100 mM Borate buffered saline, pH 8.2. No preservatives or amine-containing buffer salts added.
Concentration:	1 mg/ml
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)