

Product datasheet for **AM08206RC7-N**

Mouse IgM Isotype Control

Product data:

Product Type:	Isotype Controls
Clone Name:	11E10
Applications:	FC
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Specificity:	This antibody is specific to Lipopolysaccharide.
Buffer:	1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Conjugation:	PE-Cy7
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)