

## Product datasheet for **AM08207AC5-N**

### Mouse IgG1 Isotype Control

#### Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody is specific to T-2 Mmycotoxin.
Buffer:	1.0 mL of PBS containing 0.09% Sodium Azide as preservative and a stabilizing agent.
Conjugation:	APC-Cy5.5
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b> This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)