

Product datasheet for **AM08207C5-N**

Mouse IgG1 Isotype Control

Product data:

Product Type:	Isotype Controls
Applications:	FC
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Specificity:	This antibody is specific to T-2 Mmycotoxin.
Buffer:	1.0 ml of PBS containing 0.09% Sodium Azide as preservative
Conjugation:	Cy5
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing. Shelf life: one year from despatch.



[View online »](#)